

# DAKSH PFN 130° NAIL SYSTEM



DAKSH PFN 130°

## CARE DAKSH PFN 130° Implants



**DAKSH Proximal Femoral Nail 130° (PFN)**  
Dia 9mm to 12mm, Length 240mm



**DAKSH Proximal Femoral Nail 130° (PFN)**  
- Long Right  
Dia 9mm to 12mm, Length 320mm to 440mm




**DAKSH Proximal Femoral Nail 130° (PFN)**  
- Long Left  
Dia 9mm to 12mm, Length 320mm to 440mm



**END Cap for PFN Nail**  
Length 00mm to 15mm

## CARE DAKSH PFN 130° Instruments



**CDPFN-1900.01**

**Bone Awl**



**CDPFN-1900.02**

**Tissue Protector**



**CDPFN-1900.03**

**Guide wire for PFN locking bolts dia 2.5mm**



**CDPFN-1900.04**

**Proximal Reamer dia 13.0mm**



**CDPFN-1900.05**

**Fixation sleeve for Proximal reamer**



**CDPFN-1900.06**

**Aiming Device for Care Daksh PFN with insertion handle**



**CDPFN-1900.07**

**Nail connecting bolt cannulated**



**CDPFN-1900.08**

**Spanner**



**CDPFN-1900.09**

**Protection sleeve for 8.0mm locking bolt**



**CDPFN-1900.10**

**Protection sleeve for 6.5mm locking bolt**



**CDPFN-1900.11**

**Guide wire sleeve for 8.0mm Locking bolt**



## CARE DAKSH PFN 130° Instruments



**CDPFN-1900.12**

**Guide wire sleeve for 6.5mm  
Locking bolt**



**CDPFN-1900.13**

**Trocar for 8.0mm Locking bolt**



**CDPFN-1900.14**

**Coupling T-Handle**



**CDPFN-1900.15**

**Reamer for 8.0mm Locking  
bolt**



**CDPFN-1900.16**

**Reamer for 6.5mm Locking  
bolt**



**CDPFN-1900.17**

**Protection sleeve for 4.9mm  
locking bolt**



**CDPFN-1900.18**

**Drill sleeve for 4.9mm locking  
bolt**



**CDPFN-1900.19**

**Extractor Rod**



**CDPFN-1900.20**

**Handle For Extractor**



**CDPFN-1900.21**

**Screw driver dia 5.0mm tip for  
8.0mm & 6.5mm locking bolt**



**CDPFN-1900.22**


**Screw driver dia 3.5mm tip for  
4.9mm locking bolt**



## CARE DAKSH PFN 130° Instruments

**CDPFN-1900.23**

**Round Hammer**



**CDPFN-1900.24**

**Socket Wrench 10mm**



**CDPFN-1900.25**

**Pin Wrench**



**CDPFN-1900.26**

**Jig Assembly Connecting Parts**



**CDPFN-1900.27**

**Hammering Device**




**CDPFN-1900.28**

**Hex Cap**




**CDPFN-1900.29**

**Guide wire 2.5mm plain**




**CDPFN-1900.30**

**Guide wire 2.5mm beaded**





## CARE DAKSH PFN 130° Instruments




Specifications, size, shape, color subject to change without notice.  
Picture shape if similar to any other manufacturers is just a co-incident.

## CARE DAKSH PFN 130° Instruments

Ref. No.	Description	Material	1 SET
<b>CDPFN-1900.00</b>	<b>DAKSH PFN Nail Instrument Set</b>		
	<b>CONTENTS:-</b>		
<b>CDPFN-1900.01</b>	Bone Awl	SS	1
<b>CDPFN-1900.02</b>	Tissue Protector	SS	1
<b>CDPFN-1900.03</b>	Guide wire for PFN locking bolts dia 2.5mm	SS	1
<b>CDPFN-1900.04</b>	Proximal Reamer dia 13.0mm	SS	1
<b>CDPFN-1900.05</b>	Fixation sleeve for Proximal reamer	SS	1
<b>CDPFN-1900.06</b>	Aiming Device for Care Daksh PFN with insertion handle	SS	1
<b>CDPFN-1900.07</b>	Nail connecting bolt cannulated	SS	1
<b>CDPFN-1900.08</b>	Spanner	SS	1
<b>CDPFN-1900.09</b>	Protection sleeve for 8.0mm locking bolt	SS	1
<b>CDPFN-1900.10</b>	Protection sleeve for 6.5mm locking bolt	SS	1
<b>CDPFN-1900.11</b>	Guide wire sleeve for 8.0mm Lockinbg bolt	SS	1
<b>CDPFN-1900.12</b>	Guide wire sleeve for 6.5mm Lockinbg bolt	SS	1
<b>CDPFN-1900.13</b>	Trocar for 8.0mm Locking bolt	SS	1
<b>CDPFN-1900.14</b>	Coupling T-Handle	SS	1
<b>CDPFN-1900.15</b>	Reamer for 8.0mm Locking bolt	SS	1
<b>CDPFN-1900.16</b>	Reamer for 6.5mm Locking bolt	SS	1
<b>CDPFN-1900.17</b>	Protection sleeve for 4.9mm locking bolt	SS	1
<b>CDPFN-1900.18</b>	Drill sleeve for 4.9mm locking bolt	SS	1
<b>CDPFN-1900.19</b>	Extractor Rod	SS	1
<b>CDPFN-1900.20</b>	Handle For Extractor	SS	1
<b>CDPFN-1900.21</b>	Screw driver dia 5.0 tip for 8.0 & 6.5mm locking bolt	SS	1
<b>CDPFN-1900.22</b>	Screw driver dia 3.5mm tip for 4.9mm locking bolt	SS	1
<b>CDPFN-1900.23</b>	Round Hammer	SS	1
<b>CDPFN-1900.24</b>	Socket Wrench 10mm	SS	1
<b>CDPFN-1900.25</b>	Pin Wrench	SS	1
<b>CDPFN-1900.26</b>	Jig Assembly Connecting Parts	SS	1
<b>CDPFN-1900.27</b>	Hammering Device	SS	1
<b>CDPFN-1900.28</b>	Hex Cap	SS	1
<b>CDPFN-1900.29</b>	Guide wire 2.5mm plain	SS	1
<b>CDPFN-1900.30</b>	Guide wire 2.5mm beaded	SS	1

[www.orthocare.in](http://www.orthocare.in)  
[www.orthoimplants.com](http://www.orthoimplants.com)



FROM **WHEELS** TO **HEELS**

**ORTHO CARE**

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi – 110015, INDIA  
Phone : +91- 9810038246 | +91-11-41424511 | 41424522  
E-mail : orthocareindia2009@gmail.com | Web : [www.orthocare.in](http://www.orthocare.in)  
Factory : A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

USFDA 510K CLEARED PRODUCTS  
CE 1023 CERTIFIED PRODUCTS  
ISO 13485:2016 CERTIFICATION  
LICENSE NO. MFG/MD/2021/000347